

# bacteria spades assembly

date: 20260108

by Yanyan Liu

## what for

This is the nextflow pipeline that does bacteria genome assembly using spades. The assembly metrics are collected by QUAST. It is currently does not have the downsample option.

## how to run

### 1.copy the code below to nextflow\_spades.sh

```
#!/bin/bash

run=20251218_MiSeqi100-Morty
plates=REL_25J-B_FASTQ
sample_map=ecoli_map_run_20251218_MiSeqi100_Morty.csv
analysis=Ecoli_spades_test

/software/nextflow-align/nextflow run \
/software/nextflow-assembly-spades/main.nf \
-work-dir s3://seqwell-analysis/$run/$analysis/work \
--analysis $analysis \
--run $run \
--plates $plates \
--sample_map $sample_map \
-bg -resume
```

### 2.run as *bash nextflow\_spades.sh*

### 3.sample map example and requirement

```
sample_id,ref
REL_25J-B_A08,ecoli_REL606
REL_25J-B_B08,ecoli_REL606
REL_25J-B_C08,ecoli_REL606
REL_25J-B_D08,ecoli_REL606
REL_25J-B_E08,ecoli_REL606
REL_25J-B_F08,ecoli_REL606
REL_25J-B_G08,ecoli_REL606
REL_25J-B_H08,ecoli_REL606
```

- The sample map needs two columns, with the right header: *sample\_id* and *ref*
- For *sample\_id*, it needs to be matched to the fastq file names. For example, if the fastq file is REL\_25J-B\_A08\_R1\_001.fastq.gz and REL\_25J-B\_A08\_R2\_001.fastq.gz. The *sample\_id* will be *REL\_25J-B\_A08*